

Joint MHLW/PMDA-USP Workshop “Role of Quality in Pharmaceuticals”

Session 4 Impurities: Mutagenic impurities and more

# **Control of nitrosamine impurities in sartan drugs**

Office of Generic Drugs, PMDA  
Masahiro Uchino


## Disclaimer

Views expressed in this presentation are solely of the speaker and do not necessarily represent those of the PMDA.


# Nitrosamines

Nitrosamines are classified as probable or possible human carcinogens and referred to as “cohort of concern” compounds in the ICH M7 (R1) guideline.


## Examples of nitrosamines




*N*-Nitrosodimethylamine  
(NDMA)



*N*-Nitrosodiethylamine  
(NDEA)




*N*-Nitrosodibutylamine  
(NDBA)




*N*-Nitrosomethylphenylamine  
(NMPA)


*N*-Nitrosodiisopropylamine  
(NDIPA)



*N*-Nitrosoisopropylethylamine  
(NIPEA)



*N*-Nitroso-*N*-methyl-4-aminobutyric Acid  
(NMBA)



## Acceptable intake limits

The AI limit is a daily exposure to a nitrosamine impurity that approximates a 1:100,000 cancer risk after lifetime exposure.

The AI limits may be updated based on additional information.

Nitrosamine	AI limit (µg/day)
NDMA	0.0959
NDEA	0.0265
NDBA	0.0265
NDIPA	0.0265
NIPEA	0.0265
NMBA	0.0959

# Voluntary recalls in Japan

July 5, 2018

VALSARTAN Tablets 20 mg, 40 mg, 80 mg, 160 mg  
“AA”

- Containing valsartan API manufactured by Zhejiang Huahai Pharmaceuticals Co., Ltd.

February 7, 2019

AMVALO Combination Tablets “Pfizer”

- NDEA levels above the recommended AI limit and the detection of NDMA

September 26, 2019 -

RANITIDINE Tablets 75 mg, 150 mg (9 companies)

RANITIDINE Injection 50 mg, 100 mg (3 companies)

- The detection or possibility of the presence of NDMA

October 23, 2019

NIZATIDINE Capsules “OHARA”

- NDMA levels above the recommended AI limit

April 27, 2020

METGLUCO Tablets 250 mg, 500 mg

METFORMIN HYDROCHLORIDE Tablets 500 mg  
MT “JG”

- NDMA levels above the recommended AI limit

September 16, 2020

METFORMIN HYDROCHLORIDE Tablets 500 mg

MT “TOWA” / “NICHIIKO”

- NDMA levels above the recommended AI limit

April 26, 2021

METFORMIN HYDROCHLORIDE Tablets 500 mg

MT “JG”


- NDMA levels above the recommended AI limit or possibility of the presence of NDMA

## Root causes

Formation of nitrosamines is possible in the presence of secondary or tertiary amines and nitrites under acidic conditions.

- API processing under the presence of secondary or tertiary amines and nitrites (e.g., sartans).
- Some drug products may degrade during storage, resulting in the formation of nitrosamines (e.g., ranitidine?) .
- Excipients or container-closure system may contain amines and potential sources of nitrosating agent (e.g., metformin?) .

# Formation of nitrosamines in sartan drugs



Amine source	NOx source	Nitrosamine
<i>N,N</i> -dimethylformamide (DMF)	Sodium nitrite (NaNO <sub>2</sub> )	NDMA
<i>N</i> -methylpyrrolidone (NMP)		NMBA
Triethylamine (TEA)		NDEA
Diisopropylethylamine (DIPEA)		NDIPA
		NIPEA
Tetrabutylammonium bromide (TBAB)		NDBA

# Risk assessment in sartan drugs

## 1. Sources of secondary and tertiary amines

- If the reactions are performed at high temperatures in solvents with high boiling points (e.g., **DMF**, **NMP**) to speed up reactions, these amid solvents can degrade into secondary amines.
    - ➔ React with nitrous acid to form nitrosamines (e.g., NDMA, NMBA).
  - Tertiary amines (e.g., **TEA**, **DIPEA**) and quaternary ammonium salts (e.g., **TBAB**) used as reagents may contain secondary and tertiary amines.
    - ➔ React with nitrous acid to form nitrosamines (e.g., NDEA, NDIPA/NIPEA, NDBA).
- ✓ Nitrosamine impurities vary depending on the raw materials used in the manufacturing process.
  - ✓ In many cases, there is a risk of forming two or more nitrosamines.



# Risk assessment in sartan drugs (continued)

## 2. Presence of sodium nitrite

- Sodium nitrite is a known impurity in sodium azide, which is often used for synthesizing tetrazoles.
- An aqueous solution of sodium nitrite is used to quench residual azide, but sometimes omitted from reaction schemes and the process description in CTD module S.2.2.
- Nitrites used as reagents in one step (incl. starting materials) can carry over into subsequent steps, despite purification operations, and react with secondary or tertiary amines to generate nitrosamine impurities.

✓ The manufacturing processes of the starting materials should also be included in the risk assessment.

# Risk assessment in sartan drugs (continued)

## 3. Recovery of raw materials

- Nitrosamine impurities may be concentrated in the recovery of raw materials (e.g., solvents, reagents, and catalysts).
- Recovery of raw materials is often outsourced to third-party contractors.

### Case studies

- O-xylene and toluene were contaminated during recovery due to inadequate cleaning and to use of shared storage equipment between different customers.
- The catalyst tri-*N*-butyltin chloride was contaminated at a third-party contractor facility due to the combination of this catalyst from different customers.

# Analysis and control of nitrosamines in sartan drugs

- Manufacturers of APIs and drug products should use methods with limits of quantitation (LOQs) at or below 0.03 ppm.
- Given existing uncertainties regarding nitrosamine impurities and their presence in drugs, for at-risk APIs, testing of each batch on release should be conducted.
- If there is a risk of forming multiple nitrosamine impurities, total nitrosamine impurities should also be limited.
- Alternate approaches (e.g., upstream test of an intermediate) may be accepted based on sufficient process understanding.

Thank you for your attention.