

EXHIBIT I



调查报告 Investigation Report	编号 No.: DC _E -18003	
	版本号 Version No.: 02	页码 Page: 36 of 236

表 3-11: 缬沙坦 NDMA 检测数据汇总表

Table 3-11. Summary of NDMA test results of Valsartan

车间名称 Workshop	产品代号 Product Code	工艺 Process	标准 Spec.	检测批次 Batches	结果范围 ppm Range (ppm)	平均值 ppm Avg.(ppm)
W02 车间 Workshop W02	D5191	氯化锌 ZnCl ₂	EDMF/USDMF	1568	0.3-240.1	65.1
W02 车间 Workshop W02	D5195	三乙胺 TEA	EDMF/USDMF	64	0-7.7	0.9
2 车间 Workshop 2	C5354	三乙胺 TEA	EDMF/USDMF	122	0-37.3	1.8
2 车间 Workshop 2	C5355	氯化锌 ZnCl ₂	EDMF/USDMF	1579	0.4-136.3	63.4
2 车间 Workshop 2	C5356	三乙胺/氯化锌 TEA/ ZnCl ₂	CP	195	8.7-138.9	65.8
4 车间 Workshop 4	C5069	三乙胺 TEA	EDMF/USDMF	338	0-73.9	5.5
12/13 车间 Workshop 12/13	C5271	氯化锌 ZnCl ₂	USP	987	0-166.7	70.8
13 车间 Workshop 13	C5523	氯化锌 ZnCl ₂	CEP/USDMF	2157	0-188.1	56.7

3.1.8 其他产品基因毒性杂质 NDMA 的系统排查 Further test plan of NDMA in other products

3.1.8.1 其他产品NDMA基因毒性杂质排查计划 Test plan of NDMA in other products

- 1) 公司沙坦类产品 NDMA 基毒排查已完成，见 RARC-20180824 其他沙坦类产品潜在杂质 NDMA 的风险评估报告。

Test of NDMA in other sartans product of Huahai has been completed; refer to RARC-20180824 Quality Risk Assessment Report of NDMA in Other Sartans for details.